

## Product Datasheet

# Recombinant Mouse Acrp30/Adipoq Protein, His Tag (orb1552978)

<b>Description</b>	Recombinant Mouse Acrp30/Adipoq Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	Adipoq
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Adipoq(Glu18-Asn247) was fused with the N-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<a href="#">Q60994</a>
<b>Expiration Date</b>	6 months from date of receipt.