

Product Datasheet

Recombinant Human 4-1BBL/TNFSF9/CD137L Protein, His Tag (orb1552989)

Description	Recombinant Human 4-1BBL/TNFSF9/CD137L Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TNFSF9
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNFSF9(Arg71-Glu254) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P41273
Expiration Date	6 months from date of receipt.