

Product Datasheet

Recombinant Human Biotinylated Interleukin6(IL6) Biotinylated Protein, His Tag (orb1552994)

Description	Recombinant Human Biotinylated Interleukin6(IL6) Biotinylated Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	IL6
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human IL6(Met1-Met212) was fused with the C-terminal His Tag. The expressed protein was biotinylated in vivo by the Biotin-Protein ligase which is co-expressed and expressed in Mammalian cells.
Uniprot ID	P05231
Expiration Date	6 months from date of receipt.