

Product Datasheet

Recombinant Human EN1/Engrailed1 Protein, His tag (orb1553000)

| | |
|------------------------|--|
| Description | Recombinant Human EN1/Engrailed1 Protein, His tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | EN1 |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human EN1(Met1-Glu392) was fused with the C-terminal His Tag and expressed in E. coli. |
| Uniprot ID | Q05925 |
| Expiration Date | 6 months from date of receipt. |