

## Product Datasheet

# Recombinant Human IL12A/IL12B heterodimer Protein, Strep tag & His Tag (orb1553004)

<b>Description</b>	Recombinant Human IL12A/IL12B heterodimer Protein, Strep tag & His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	IL12A/IL12B
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IL12A(Met1-Ser219) was fused with the C-terminal Strep Tag. This protein contains the human IL12B(Met1-Ser328) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">P29460</a> , <a href="#">P29459</a>
<b>Expiration Date</b>	6 months from date of receipt.