

Product Datasheet

Recombinant General Histone H3 Protein, His Tag (orb1553005)

Description	Recombinant General Histone H3 Protein, His Tag
Species/Host	E. coli
Reactivity	All
Endotoxins	Please contact our customer service for more information.
Target	HIST1H3A
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the general HIST1H3A(Met1-Ala136)with C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q6LBF0 , Q6LED0 , P68431 , P68433
Expiration Date	6 months from date of receipt.