

Product Datasheet

Recombinant *Shigella flexneri* Thioredoxin/trxA Protein, His Tag (orb1553007)

Description	Recombinant <i>Shigella flexneri</i> Thioredoxin/trxA Protein, His Tag
Species/Host	<i>E. coli</i>
Reactivity	Bacteria
Endotoxins	Please contact our customer service for more information.
Target	trxA
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the <i>Shigella flexneri</i> trxA(Met1-Ala109) was fused the C-terminal His Tag and expressed in <i>E. coli</i> .
Uniprot ID	P0AA30
Expiration Date	6 months from date of receipt.