

## Product Datasheet

# Recombinant Mouse IL12b Protein, TwinStrep tag (orb1553011)

<b>Description</b>	Recombinant Mouse IL12b Protein, TwinStrep tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	IL12b
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse IL12b(Met1-Ser335) was fused with the C-terminal Twin-Strep Tag and expressed in Mammalian.
<b>Uniprot ID</b>	<a href="#">P43432</a>
<b>Expiration Date</b>	6 months from date of receipt.