

Product Datasheet

Recombinant Mouse Il12a&Il12b heterodimer Protein, His Tag & TwinStrep tag (orb1553012)

Description	Recombinant Mouse Il12a&Il12b heterodimer Protein, His Tag & TwinStrep tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Il12a/Il12b
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	<p>This protein contains the mouse Il12a(Met1-Ala215) was fused with the C-terminal His Tag.</p> <p>This protein contains the mouse Il12b(Met1-Ser335) was fused with the C-terminal Twin-Strep Tag and expressed in Mammalian.</p>
Uniprot ID	P43431, P43432
Expiration Date	6 months from date of receipt.