

Product Datasheet

Recombinant Human TPO Protein, His Tag (orb1553019)

Description	Recombinant Human TPO Protein, His Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TPO
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human Tpo(Met1-Arg846) was fused with the C-terminal His Tag and expressed in Baculovirus-Insect cells.
Uniprot ID	P07202
Expiration Date	6 months from date of receipt.