

Product Datasheet

Recombinant Mouse Padi2/Pad2 Protein, TwinStrep Tag (orb1553020)

Description	Recombinant Mouse Padi2/Pad2 Protein, TwinStrep Tag
Species/Host	Virus
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Padi2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Padi2(Met1-Pro673) was fused with the N-terminal TwinStrep Tag and expressed in Baculovirus-Insect cells.
Uniprot ID	Q08642
Expiration Date	6 months from date of receipt.