

Product Datasheet

Recombinant MoMLV gag-pol Protein, His Tag (orb1553021)

Description	Recombinant MoMLV gag-pol Protein, His Tag
Species/Host	E. coli
Reactivity	Virus
Endotoxins	Please contact our customer service for more information.
Target	gag-pol
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the moloney murine leukemia virus gag-pol(Thr660-Leu1330) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P03355
Expiration Date	6 months from date of receipt.