

## Product Datasheet

# Recombinant Human PADI2 Protein, TwinStrep Tag (orb1553035)

<b>Description</b>	Recombinant Human PADI2 Protein, TwinStrep Tag
<b>Species/Host</b>	Virus
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	PADI2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human PADI2(Met1-Pro665) was fused with the N-terminal TwinStrep Tag and expressed in Baculovirus-Insect cells.
<b>Uniprot ID</b>	<a href="#">Q9Y2J8</a>
<b>Expiration Date</b>	6 months from date of receipt.