

Product Datasheet

Recombinant Human PDGFB Protein, His Tag (orb1553037)

Description	Recombinant Human PDGFB Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	PDGFB
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PDGFB(Glu21-Ala241) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P01127
Expiration Date	6 months from date of receipt.