

Product Datasheet

Recombinant Human DHX36 Protein, Gst Tag (orb1553056)

Description	Recombinant Human DHX36 Protein, Gst Tag
Species/Host	Virus
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	DHX36
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human DHX36(Ser2-Ser1008) was fused with the N-terminal Gst Tag and expressed in Baculovirus-Insect cells
Uniprot ID	Q9H2U1
Expiration Date	6 months from date of receipt.