

Product Datasheet

Recombinant Human TNNC1&TNNI3 Protein, His Tag (orb1553071)

Description	Recombinant Human TNNC1&TNNI3 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	TNNC1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNNI3(Met1-Glu161 & Thr31-Arg111) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P19429, P63316
Expiration Date	6 months from date of receipt.