

Product Datasheet

Recombinant Human CD40L/CD40LG Protein, His Tag (orb1553087)

Description	Recombinant Human CD40L/CD40LG Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CD40LG
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CD40LG(Glu108-Leu261) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P63304
Expiration Date	6 months from date of receipt.