

Product Datasheet

Recombinant Human GITRL/TNFSF18 Protein, His Tag (orb1553102)

| | |
|------------------------|---|
| Description | Recombinant Human GITRL/TNFSF18 Protein, His Tag |
| Species/Host | Virus |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | TNFSF18 |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human TNFSF18(Glu74-Ser199) was fused with the N-terminal His Tag and expressed in Baculovirus-Insect cells |
| Uniprot ID | Q9UNG2 |
| Expiration Date | 6 months from date of receipt. |