

## Product Datasheet

# Recombinant Human FOLH1 Protein, TwinStrep Tag (orb1553110)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human FOLH1 Protein, TwinStrep Tag  |
| <b>Species/Host</b>    | Virus   |
| <b>Reactivity</b>      | Human   |
| <b>Endotoxins</b>      | Please contact our customer service for more information.   |
| <b>Target</b>          | FOLH1   |
| <b>Note</b>            | For research use only   |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE  |
| <b>Source</b>          | This protein contains the human FOLH1(Lys44-Ala750) was fused with the N-terminal TwinStrep Tag and expressed in Baculovirus-Insect cells |
| <b>Uniprot ID</b>      | <a href="#">Q04609</a>  |
| <b>Expiration Date</b> | 6 months from date of receipt.  |