

Product Datasheet

Recombinant Human CD334/FGFR4 Protein, His Tag (orb1553147)

Description	Recombinant Human CD334/FGFR4 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	FGFR4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human FGFR4(Met1-Asp369) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P22455
Expiration Date	6 months from date of receipt.