

Product Datasheet

Recombinant Human CD333/FGFR3 Protein, His Tag (orb1553148)

| | |
|------------------------|--|
| Description | Recombinant Human CD333/FGFR3 Protein, His Tag |
| Species/Host | Mammal |
| Reactivity | Human |
| Endotoxins | < 1.0 EU per µg of the protein as determined by the LAL method. |
| Target | FGFR3 |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human FGFR3(Met1-Gly377) was fused with the C-terminal His Tag and expressed in Mammalian cells. |
| Uniprot ID | P22607 |
| Expiration Date | 6 months from date of receipt. |