

Product Datasheet

Recombinant Human CD328/SIGLEC7 Protein, Fc Tag (orb1553152)

| | |
|------------------------|--|
| Description | Recombinant Human CD328/SIGLEC7 Protein, Fc Tag |
| Species/Host | Mammal |
| Reactivity | Human |
| Endotoxins | < 1.0 EU per µg of the protein as determined by the LAL method. |
| Target | SIGLEC7 |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human SIGLEC7(Gln19-Gly357) was fused with the C-terminal Fc Tag and expressed in Mammalian cells. |
| Uniprot ID | Q9Y286 |
| Expiration Date | 6 months from date of receipt. |