

## Product Datasheet

# Recombinant Human CD254/TNFSF11 Protein, His Tag (orb1553216)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human CD254/TNFSF11 Protein, His Tag  |
| <b>Species/Host</b>    | Mammal  |
| <b>Reactivity</b>      | Human   |
| <b>Endotoxins</b>      | < 1.0 EU per µg of the protein as determined by the LAL method.   |
| <b>Target</b>          | TNFSF11   |
| <b>Note</b>            | For research use only   |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE  |
| <b>Source</b>          | This protein contains the human TNFSF11(Tyr69-Asp317) was fused with the C-terminal His Tag and expressed in Mammalian cells. |
| <b>Uniprot ID</b>      | <a href="#">O14788</a>  |
| <b>Expiration Date</b> | 6 months from date of receipt.  |