

Product Datasheet

Recombinant Human CD218B/IL18RAP Protein, His Tag (orb1553241)

Description	Recombinant Human CD218B/IL18RAP Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	IL18RAP
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human IL18RAP(Met1-Gly357) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	O95256
Expiration Date	6 months from date of receipt.