

Product Datasheet

Recombinant Human CD217/IL17RA Protein, His Tag (orb1553243)

| | |
|------------------------|--|
| Description | Recombinant Human CD217/IL17RA Protein, His Tag |
| Species/Host | Mammal |
| Reactivity | Human |
| Endotoxins | < 1.0 EU per µg of the protein as determined by the LAL method. |
| Target | IL17RA |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human IL17RA(Met1-Trp 320) was fused with the C-terminal His Tag and expressed in Mammalian cells. |
| Uniprot ID | Q96F46 |
| Expiration Date | 6 months from date of receipt. |