

## Product Datasheet

# Recombinant Human CD159A/KLRC1 Protein, His Tag (orb1553280)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Human CD159A/KLRC1 Protein, His Tag  |
| <b>Species/Host</b>    | Mammal   |
| <b>Reactivity</b>      | Human  |
| <b>Endotoxins</b>      | < 1.0 EU per µg of the protein as determined by the LAL method.  |
| <b>Target</b>          | KLRC1  |
| <b>Note</b>            | For research use only  |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE   |
| <b>Source</b>          | This protein contains the human KLRC1(Arg100-Leu233) was fused with the N-terminal His Tag and expressed in Mammalian cells. |
| <b>Uniprot ID</b>      | <b>P26715</b>  |
| <b>Expiration Date</b> | 6 months from date of receipt.   |