

Product Datasheet

Recombinant Human CD140B/PDGFRB Protein, His Tag (orb1553306)

| | |
|------------------------|---|
| Description | Recombinant Human CD140B/PDGFRB Protein, His Tag |
| Species/Host | Mammal |
| Reactivity | Human |
| Endotoxins | < 1.0 EU per µg of the protein as determined by the LAL method. |
| Target | PDGFRB |
| Note | For research use only |
| Purity | ≥90% as determined by SDS-PAGE |
| Source | This protein contains the human PDGFRB(Met1-Lys531) was fused with the C-terminal His Tag and expressed in Mammalian cells. |
| Uniprot ID | P09619 |
| Expiration Date | 6 months from date of receipt. |