

## Product Datasheet

# Recombinant Human CD120A/TNFRSF1A Protein, His Tag (orb1553328)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Human<br>CD120A/TNFRSF1A Protein, His Tag  |
| <b>Species/Host</b>    | Mammal   |
| <b>Reactivity</b>      | Human  |
| <b>Endotoxins</b>      | < 1.0 EU per µg of the protein as<br>determined by the LAL method.   |
| <b>Target</b>          | TNFRSF1A   |
| <b>Note</b>            | For research use only  |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE   |
| <b>Source</b>          | This protein contains the human<br>TNFRSF1A(Met1-Thr211) was fused<br>with the C-terminal His Tag and<br>expressed in Mammalian cells. |
| <b>Uniprot ID</b>      | <b>P19438</b>  |
| <b>Expiration Date</b> | 6 months from date of receipt.   |