

## Product Datasheet

# Recombinant Human CD112/NECTIN2 Protein, His Tag (orb1553336)

<b>Description</b>	Recombinant Human CD112/NECTIN2 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	NECTIN2
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human NECTIN2(Met1-Gly360) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q92692">Q92692</a>
<b>Expiration Date</b>	6 months from date of receipt.