

Product Datasheet

Recombinant Human CD85A/LILRB3 Protein, His Tag (orb1553373)

Description	Recombinant Human CD85A/LILRB3 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	LILRB3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human LILRB3(Met1-Glu443) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	O75022
Expiration Date	6 months from date of receipt.