

## Product Datasheet

# Recombinant Human CD49A/ITGA1 Protein, His Tag (orb1553416)

|                        |                                                                                                                             |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human CD49A/ITGA1 Protein, His Tag                                                                              |
| <b>Species/Host</b>    | Mammal                                                                                                                      |
| <b>Reactivity</b>      | Human                                                                                                                       |
| <b>Endotoxins</b>      | < 1.0 EU per µg of the protein as determined by the LAL method.                                                             |
| <b>Target</b>          | ITGA1                                                                                                                       |
| <b>Note</b>            | For research use only                                                                                                       |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE                                                                                              |
| <b>Source</b>          | This protein contains the human ITGA1(Met1-Pro1141) was fused with the C-terminal His Tag and expressed in Mammalian cells. |
| <b>Uniprot ID</b>      | <b>P56199</b>                                                                                                               |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                              |