

## Product Datasheet

# Recombinant Human IL23A/IL-23 Protein, His Tag (orb1553488)

<b>Description</b>	Recombinant Human IL23A/IL-23 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	IL23A
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IL23A(Met1-Pro189) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">Q9NPF7</a>
<b>Expiration Date</b>	6 months from date of receipt.