

Product Datasheet

Recombinant Human IL6 Protein, His tag (orb1553506)

Description	Recombinant Human IL6 Protein, His tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	IL6
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human IL6(Met1-Met212) was fused with the C-terminal His tag and expressed in Mammalian cells.
Uniprot ID	P05231
Expiration Date	6 months from date of receipt.