
Product Datasheet

Acetyl CoA Carboxylase(Phospho-Ser79) antibody (orb155573)

Description

Acetyl CoA Carboxylase(Phospho-

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from mouse Acetyl Coenzyme A Carboxylase around the phosphorylation site of Ser79 (SM(p-S)GL)
Target	Phospho-ACACA (Ser79)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	266 kDa
Uniprot ID	Q13085
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.