

Product Datasheet

alpha Adducin (Phospho-Thr445) antibody (orb155638)

Description	alpha Adducin (Phospho-Thr445) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human alpha Adducin around the phosphorylation site of Thr445 (EK(p-T)RW)
Target	phospho-alpha Adducin (Thr445)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80 kDa
Uniprot ID	P35611
Expiration Date	12 months from date of receipt.