



Product Datasheet

APC (Phospho-Ser2054) antibody (orb155681)

Description	APC (Phospho-Ser2054) antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human APC around the phosphorylation site of Ser2054 (KP(p-S)RL)
Target	phospho-APC (Ser2054)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	312 kDa
Uniprot ID	P25054
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse colon tissue using p...



IHC-P analysis of rat cerebellum using p...