



Product Datasheet

APOBEC3B antibody (orb155694)


Description	APOBEC3B antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human APOBEC3B. Please contact us for the exact immunogen sequence. The peptide is available as orb12939.
Target	APOBEC3B
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	-100 µg (in 200 µl): 0.5 mg/ml-200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	46 kDa
Uniprot ID	F7EWS7 , Q9UHL7
NCBI	001270411 , 001257340
Entrez	9582
Dilution Range	IF/ICC: 1:100-600, IHC-P: 1:100-600, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



IF image of rat testis tissue using APOB...



IHC-P image of rat testis tissue using A...



Western blot analysis of rat lymph node ...