

Product Datasheet

ATRIP (Phospho-Ser68+Ser72) antibody (orb155766)

Description	ATRIP (Phospho-Ser68+Ser72) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human ATRIP around the phosphorylation site of Ser68 + Ser72 (LA(p-S)QAL(p-S)QC)
Target	phospho-ATRIP (Ser68+Ser72)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86 kDa
Uniprot ID	Q8WXE1
Expiration Date	12 months from date of receipt.