



Product Datasheet

c-ABL(Phospho-Tyr115) antibody (orb156230)

Description	c-ABL(Phospho-Tyr115) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human c-Abl (SN(p-Y)IT)
Target	phospho-c-ABL (Tyr115)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	P00519
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of rat brain using phospho-c-ABL (Tyr115) antibody.



Western blot: Lane 1: SiHa (Human) Cell; Lane 2: La...