


## Product Datasheet

# Flag Tag (3A6) Mouse mAb (HRP) Antibody (orb1562806)

<b>Description</b>	Flag Tag (3A6) Mouse mAb (HRP) Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Other
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	WB
<b>Immunogen</b>	DYKDDDDK synthetic peptide conjugated to KLH.
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	WB: 1/1000-1/10000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Flag Tag (3A6) ...