


Product Datasheet

GST Tag Rabbit pAb (HRP) Antibody (orb1562807)

Description	GST Tag Rabbit pAb (HRP) Antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	HRP
Tested Applications	WB
Immunogen	Purified recombinant Escherichia coli BL21(DE3) GST protein expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 26 kDa; Observed MW:
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of purified recomb...