

Product Datasheet

His Tag Rabbit pAb (HRP) Antibody (orb1562808)

Description	His Tag Rabbit pAb (HRP) Antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	HRP
Tested Applications	WB
Immunogen	His-tag Rabbit pAb (HRP-Conjugate) is produced by immunizing animals with a 6xHis synthetic peptide.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.