


Product Datasheet

Dsba Antibody (orb1562843)

Description	Dsba Antibody
Species/Host	Mouse
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant E.coli. DsbA protein fragments over expressed in E.coli.
Target	Dsba
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	4.3 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	Calculated MW: 23 kDa; Observed MW:
Uniprot ID	POAEG4
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of DsbA with Dsba ...