

Product Datasheet

c-Myc Tag Rabbit pAb Antibody (orb1562846)

Description	c-Myc Tag Rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide EQKLISEEDL coupled to KLH.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB: 1/1000-1/10000
Expiration Date	12 months from date of receipt.



1ug c-Myc
fusion
protein+
Primary
antibo...