


Product Datasheet

Flag Tag Rabbit mAb Antibody (orb1562848)

Description	Flag Tag Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide of DDDDK tag
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Form/Appearance	Liquid
Concentration	0.3 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB: 1/1000-1/10000 IHC: 1/50-1/200 IF: 1/100-1/500 IP: 1/20
Expiration Date	12 months from date of receipt.



Western blot analysis of Flag Tag in CHO...