



## Product Datasheet

### TUBA1B Antibody (orb1562970)


<b>Description</b>	TUBA1B Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC-Fr, IHC-P, IP, WB
<b>Immunogen</b>	Recombinant Protein of Tubulin alpha-1A chain
<b>Target</b>	TUBA1B
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	3.4 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	Calculated MW: 50 kDa; Observed MW: 50 kDa
<b>Uniprot ID</b>	<b>P68363</b>
<b>Dilution Range</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of alpha Tubulin i...



Immunofluorescence analysis of alpha Tub...



Immunoprecipitation analysis of alpha Tu...