



## Product Datasheet

### ENO2 Antibody (orb1563158)


<b>Description</b>	ENO2 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide derived from human NSE.
<b>Target</b>	ENO2
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	4.5 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	Calculated MW: 47 kDa; Observed MW: 47 kDa
<b>Uniprot ID</b>	<b>P09104</b>
<b>Dilution Range</b>	WB: 1/500-1/1000 ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of NSE (1H8) in ra...



Western blot analysis of Enolase2 in SH-...



Observed Enolase2 levels in CHO-K1 cell

...