



Product Datasheet

JAK1 Antibody (orb1563180)


Description	JAK1 Antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Purified recombinant human Jak1 protein fragments expressed in E.coli.
Target	JAK1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	1.7 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	Calculated MW: 133 kDa; Observed MW: 130 kDa
Uniprot ID	P23458
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20
Expiration Date	12 months from date of receipt.



Western blot analysis of JAK1 (8B8) in C...



Immunocytochemistry analysis of JAK1 (8B...



Immunoprecipitation analysis of JAK1 (8B...