



Product Datasheet

EIF4EBP1 Antibody (orb1563205)


Description	EIF4EBP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide.
Target	EIF4EBP1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	0.928 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Calculated MW: 13 kDa; Observed MW: 15-20 kDa
Uniprot ID	Q13541
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.



Immunoprecipitation analysis of 4E BP1 i...



Western blot analysis of 4E BP1 in SH-SY5...



Western blot analysis of 4E BP1 in Jurka...