



Product Datasheet

BRAP Antibody (orb1563320)

Description	BRAP Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human BRAP protein fragments expressed in E.coli.
Target	BRAP
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Form/Appearance	Liquid
Concentration	4 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	Calculated MW: 67 kDa; Observed MW: 67 kDa
Uniprot ID	Q7Z569
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.



Western blot analysis of BRAP (1E7) in C...



Western blot analysis of BRAP (1E7) in A...